



CAMEL HOPE-YE246

# BRUSHED WHEELCHAIR

This Power wheelchair is **made of Aircraft level aluminum material**. It is our first innovation which putting Brush motor on such a lightweight and simple designed power wheelchair. Huge improvements made on the innovated structure, appearance and product quality. Mean while, the competitive price let more passengers affordable for this unique wheelchair.

Folding in **1s**



YATTLL INDUSTRY CO., LTD.

Tel: +86 755 817 071 99 Web: [www.yattll.com](http://www.yattll.com) Email: [info@yattll.com](mailto:info@yattll.com)

Factory: No.2 Lingchang Road, PingShan Village, TangXia Town, DongGuang City, GuangDong, China

Office: RM.1403, BLK.B of Duocai Office Complex, GuanlanHigh-tech Zone, Longhua DIST., Shenzhen, China




# FOCUS ON DETAILS



## Anti-tip

Uphill to prevent sliding or displacement after being hit by others.



## DOUBLE BAR

Double Base Bar To Ensure The Wheelchair Be Straight, And Wheels Never Go Skew.



## PU SOLID TYRE

The rear wheels have strong pattern for good grip on all surfaces.



## FOLDABLE LEGREST

The comfortable legrest can be folded up for easy entry and exit.



## ADJUSTABLE ARMREST

The two sides can be lifted and retracted freely.

**contact us**

YATTLL INDUSTRY CO., LTD.

Tel: +86 755 817 071 99 Web: [www.yattll.com](http://www.yattll.com) Email: [info@yattll.com](mailto:info@yattll.com)

Factory: No.2 Lingchang Road, PingShan Village, TangXia Town, DongGuang City, GuangDong, China

Office: RM.1403, BLK.B of Duocai Office Complex, Guanlan High-tech Zone, Longhua DIST., Shenzhen, China



# FOLDING 1S

## Parameter



**Front wheel:** 8", PU Solid Tyre

**Rear Wheel:** 12.5", PU Solid Tyre

**Net Weight:** 26.5 kg

**Brake:** Electromagnetic Braking system

**Batteries:** Lithium Battery 24V 11Ah

**Motors (pair):** Brushless Motor 250W

**Charger:** High Power 24V 2Ah

**Controller:** Reliable Brushless Controller

**Maximum Range:** 10-13 km

**Suspension System:** Yes

**Ground Clearance:** 8 cm

**Max Speed:** 8 km/h

**Charging Hour:** 5-6 Hours

**Max Slope:** 12 degrees

**Weight Capacity:** 135 KG

**Dimension (Folded):** 80\*62\*34cm

## Features


- Top seller model with classic design
- 24V 350W powerful brush motor
- Elevating Armrest, Length adjustable
- Foldable backrest, detachable footrest
- Folding without battery remove

## Footnotes

1) Range and speed vary with userweight-terrain type, battery charge battery condition and tyre pressure.

2) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.

3) Battery weight may vary based on manufacturer.



YATTLI INDUSTRY CO., LTD.

Tel: +86 755 817 071 99 Web: [www.yattli.com](http://www.yattli.com) Email: [info@yattli.com](mailto:info@yattli.com)

Factory: No.2 Lingchang Road, PingShan Village, TangXia Town, DongGuang City, GuangDong, China

Office: RM.1403, BLK.B of Duocai Office Complex, GuanlanHigh-tech Zone, Longhua DIST., Shenzhen, China

