

## Product Parameter

Front wheel	8" wheel
Rear Wheel	12.5" wheel
Net Weight	27kg
Brake	Electromagnetic Braking system
Batteries	Lithium Battery 24V 6Ah*2
Motors (pair)	Brushless Motor 180W
Charger	24V 2A/3A/5A
Controller	Reliable Brushless Controller
Maximum Range	15-18 km
Suspension System	Yes
Ground Clearance	8 cm
Max Speed	8 km/h
Charging Hour	5-6 Hours
Max Slope	12 degrees
Weight Capacity	150 KG
Dimension (Folded)	79*60*37cm
Dimension (Folding)	110*60*98cm

1. Range and speed vary with user weight terrain type, battery charge battery condition and tyre pressure.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
3. Battery weight may vary based on manufacturer.

# TO BE A YATLL AGENT

## 2021

### Requirements

1. At least 1 lawful and registered local shop or sale place managed by yourself.
2. Legally registered company.
3. Have some knowledge and strong interest in electric wheelchairs.

### Our Support

1. Agency area protection.
2. The best quality and competitive price.
3. Prior rights to the new products with new opportunity.
4. No minimum ordering quantity.
5. Soonest response on any of your questions.
6. Free parts supply support policy.

JOIN US



**YATLL INDUSTRY CO., LTD.**



**Tel:** +86 755 8170 7199

**Web:** www.yatll.com

**Email:** info@yatll.com

**Factory:** No.1 Lingchang Road, PingShan Village, TangXia Town, DongGuang City, GuangDong, China.

**Office:** RM.1403, BLK.B of Duocai Office Complex, Gu-anlan High-tech Zone, Longhua DIST., Shenzhen, China.

New product



# Quick release Battery

Model: **YE200**



www.yatll.com |

## Detail display chart



## NEW PRODUCT

YE200 has a unique design quick release battery, It is removable and much convenient to charge. The whole wheelchair is precision machining made, the beautiful finish of the frame makes it be a high end product.

### Features

1. Quick release battery design, easy removable battery.
2. Two SAMSUNG battery available to make a longer driving distance.
3. Smart control system, soft and smooth.
4. Reinforced frame making the wheelchair maximum loading 150kg.
5. High power motor pair, silent and drives longer.
6. Ergonomically designed cushion, comfortable and removable.

### Colors

Cushion color



Frame color



### Folding steps



Three steps quick folding , Only takes 1s.

• YATLL INDUSTRY CO., LTD.

**Tel:** +86 755 8170 7199

**Email:** info@yatll.com